

NEW IMPROVED MODEL SE7500

Ideal for the slightly to moderately hearing impaired person, the SuperEarPlus® is perfect for watching TV or getting the most from the theatre, meetings, lectures, worship services or simple conversations. As an Assistive Listening Device, the SuperEarPlus is also ADA compliant.

Newly re-designed microphone and carrying case help eliminate annoying case noise caused by handling.

Unobtrusive earbuds are easily camouflaged in public. Comfortable stereo headphones (also included) are great for home use and provide hearing aid wearers with a welcome relief from in-ear aids.

Laboratory tested to hearing aid standards, the SuperEarPlus increases ambient sound gain by more than 40 decibels and picks up sounds from up to 100 yards away using a highly sensitive microphone.

Stylish design weighs just 3 ounces, measures 3.5" x 2" x .75" and fits easily into pocket or purse.

SuperEar^{Plus}® PERSONAL SOUND AMPLIFIER

Features include:

- *Laboratory Tested to Hearing Aid Standards.*
- *Volume control wheel, located on the headphone lead wires, permits easy volume adjustment.*
- *Convenient carrying case clips to pocket or belt and helps minimize case noise.*
- *45 minute Automatic Shut-off to conserve battery life. 2 AAA alkaline batteries (included) provide up to 80 hours of operation.*
- *Lightweight and unobtrusive stereo "earbuds" and comfortable headphones are both included.*
- *Lighted "On/Off" button.*
- *100 to 14000 hertz frequency range.*
- *5 Year Limited Warranty.*



SonicTechnology®

120 Richardson St., Grass Valley, CA 95945
800-247-5548 ■ (530) 272-4607 ■ www.sonictechnology.com